Electricité

Electrostatique

Electricity

Electrostatics

Réf :
272 015

Français – p 1

English – p 5

Version : 0201

Coulombmètre Initio®

Initio® coulometer

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2. Examples of experiments

❖ *Electrification by contact:*

1) Discharge the plate by pressing the Reset button.

Rub an ebonite stick with cat fur. The stick gets negatively charged. Take the end of the stick close to the charge collector plate until it is in contact with the latter.

A part of the negative charge is transferred. The plate is negatively charged. The (-) sign and the value of the charge in nC remain displayed even after the stick is removed.

2) Discharge the plate by pressing the Reset button.

Rub an ebonite stick with cat fur. The stick gets negatively charged.

Touch the end of the stick and then the entire stick to the plate.

Compare the 2 values.

(To predict the sign of the charge, refer to the triboelectricity list on the following page)

❖ *Electrification by influence:*

1) Discharge the plate by pressing the Reset button.

Rub an ebonite stick with cat fur. Take the stick close to the charge collector plate without touching it. The (-) sign and the value of the charge in nC are displayed only when the stick is close to the plate.

2) Discharge the plate by pressing the Reset button.

Rub an ebonite stick with cat fur. Take the stick close to the collector plate and press the reset button. The free electrons of the plate are repelled by negative charges of the stick and discharged by the ground. By removing the stick and the finger from the Reset button, the plate remains positively charged since it has lost electrons. The value of the charge in nC are displayed (the (+) sign is implicit).

Triboelectricity list

	Acetate
	Rabbit fur
	Glass
	Dry human hair
	Nylon
MORE POSITIVE (+)	Wool
	Cat fur
	Lead
	Silk
	Aluminium
	Paper
	Cotton
	Wood
	Amber
	Steel
	Ebonite
	Nickel - Copper
MORE NEGATIVE (-)	Brass – Stainless steel
	Rubber
	Most thermoplastics (PVC, etc.)
	PTFE (polytetrafluoroethylene, teflon)

When we rub a material **A** with material **B** placed after it in the list, **A** gets charged with positive electrical charges and **B** with negative electrical charges.

For example, glass rubbed with wool, cat's fur, cotton or nylon will always be positive. It will be negative if it is rubbed with rabbit fur.

At the top of the list are materials that have a tendency of losing electrons and at the bottom of the list are materials that have maximum "affinity" for electrons.


The intensity of the triboelectric charge depends on the following factors:

- Separation between the two materials in the triboelectric table
- Experimental conditions (temperature, relative humidity)
- Conditions of the surface of the materials (cleanliness, texture)
- Contact pressure during rubbing
- Velocity of separating the materials after rubbing.

3. Technical specifications

- Range : - 199.9 to + 199.9 nC
- Resolution : 0.1 nC
- Internal capacitance : 1 μ F – 250 V
- Dimensions : 100 x 100 x 40 mm
- Power supply : 1 6F22 type 9V battery
- Battery life : 60 hours

4. Replacing the battery

To open the battery housing, insert a small screw driver in the notch at the edge of the flap and lever it up as indicated in the drawing hereof.

5. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.